



Agile

The Agile is simple, but never ordinary. It is built with aircraft grade aluminum, which makes it light in weight and durable.





Easy Transfer

The flip-back armrests make moving from side to side easy and quick.

The 80° removable footrests ensure a safe transfer.



Smooth and Agile Movement

Agile wheelchair is designed with the carer in mind.

With advanced engineering, this chair is very easy to push and manoeuvre.



Compact and Lightweight

The Agile weighs from 12.7 KG complete, which is lighter than your carry on luggage.

With the backrest and seat folded, the chair's compact size makes it very convenient to carry and travel with.

Lightweight and Transportable

The Agile is simple, but never ordinary. Agile is built with aircraft grade aluminum, which makes it lighter in weight (12.7KG) and durable in use.



The Curved Line Design helps guide the seat upholstery into a nice and neat fold. Quick-Release Rear Wheels allow lighter carrying weight and smaller storage space.

▲ Black Frame, 18" Seat Depth, 24" Rear Wheel
With Folding Backrest

The Agile is available in a range of configurations, that will fit the body dimensions of most users.

Seat Height
505mm

Seat Width
16" 、18" 、20"

Seat Depth
410mm/ 460mm

Backrest Frame Types
Foldable

Rear Wheel Types
Fixed/ Quick-released




Quick to Adapt. Ready for Use



Attendant-Propelled Model



Self-Propelled Model




Adjustable Rear Wheel Axles

By adapting the optional rear-wheel extension plate, you can configure the position of rear wheel axles and optimise Agile's maneuverability or stability.

Robust and Reliable

The Agile is approved for transportation in accordance with ISO 7176-19 Crash Test. This means it is suitable to be used as a seat in a motor vehicle for transportation.





Optional Accessories

Tension Adjustable Backrest

Unlike traditional backrest designs, the tension adjustable backrest is made up of high-quality Velcro straps. By changing the tension of each strap, the ideal curvature can be created to fit the user's back.



Attendant Brakes



Flip-back,
Height Adjustable
Armrests




Elevating
Legrests



Anti-Tips



Seat Belt



Crutch Holder



I.V. Pole Holder



Oxygen Bottle
Holder



Modern Look

When it comes to the style, you can mix and match the frame color (black, white, and pearl silver) with the classic black or the grayish dot upholstery.



Scan to learn more about Agile

Agile Specifications

| Seat Width* Seat Depth (inch) | 16 / 18 / 20 *16 / 18 | 16 / 18 / 20 *16 / 18 | 16 / 18 / 20 *16 / 18 | 16 / 18 / 20 *16 / 18 |
|----------------------------------|--|---|---|---|
| Seat Width (mm) | 430 / 480 / 530 | 430 / 480 / 535 | 430 / 480 / 530 | 430 / 480 / 530 |
| Seat Depth (mm) | 510 / 560 | 505 / 555 | 505 / 555 | 505 / 555 |
| F/R wheels (inch) | 7 / 14 | 7 / 22 | 8 / 14 | 8 / Q24 |
| Seat Height (mm) | 480 | 480 | 505 | 505 |
| Backrest Height (mm) | 430 / 510 | 430 / 510 | 430 / 510 | 430 / 510 |
| Armrest Height (mm) | 220 | 220 | 220 | 220 |
| Calf Length (mm) | 380-480 | 380-480 | 380-480 | 380-480 |
| Max. User Weight (kg) | 125 | 125 | 125 | 125 |
| Total Weight (kg) | 12.3-12.8 12.6-13.1* 12.8-13.3** | 12.9-13.4 13.2-13.7* 13.4-13.9** | 12.2-12.7 12.5-13.0* 12.7-13.2** | 13.9-14.4 14.2-14.7* 14.4-14.9** |
| Overall L x W x H (mm) | Min. 1025*600*915 Max. 1095*700*995 | Min. 1055*610*920 Max. 1125*710*1000 | Min. 1025*600*935 Max. 1095*700*1015 | Min. 1165*605*940 Max. 1120*705*1020 |
| Overall Folded L x W x H (mm) | Min. 710*280*735 Max. 830*280*915 | Min. 710*280*765 Max. 820*280*920 | Min. 750*280*760 Max. 840*280*930 | Min. 850*280*750 Max. 1070*280*750 |

* The seat depth option for 20" seat width is only 18" .

* 430 mm foldable backrest/ ** 510 mm foldable backrest.

* The actual product specifications may vary according to different configurations.

* Karma reserves the right to modify information herein without further notice.

Unit:cm ± 1cm ± 0.5kg

Rev.02 Nov.2021



KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

F +886-5-206 7788

Emailglobalsales@karma.com.tw

soma

Designed and Made by Karma