

Agile

Designed with users and attendants in mind, Agile is the world's new standard for everyday wheelchairs. Detachable, swingaway footrests and flip-up armrests make transfers quick and easy, while the foldable frame is lightweight and compact to maximize portability and ease of use. Various configurations, options, and accessories allow Agile to fulfill the needs of a wide range of users.



Rear Wheel Extension Plate (Optional)

Adjust the horizontal position of the rear wheel to improve propelling efficiency or stability.



Tension Adjustable Backrest (Optional)

The internal structure of the backrest is formed by a series of individually-adjustable Velcro straps, which can create a user-specific fit to maximize comfort and support.



FRAME COLORS



Scan to learn more about Agile



Standard



Easy Transfer with Agile

The standard flip-back armrests make moving from side to side easy and quick. The removable 80° footrests enable safe transfers into and out of the wheelchair.



Optional



Upgrades

Optional upgrades include the flip-up, height adjustable armrests and elevating legrests.



Lightweight, easy to fold, and compact, Agile is a practical option for users and attendants, alike.

Series	Product Name	Nominal seat width (mm/in)	Nominal seat depth (mm/in)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)							
SOMA	AGILE	16	16	7 x 14	480	430 / 510	220	380 - 480	14	8							
		18															
		16	18								7 x F22	480	430 / 510	220	380 - 480	14	8
		18															
		20															
		16	16	8 x 14	505	430 / 510	220	380 - 480	14	8							
		18															
		16	18								8 x Q24	505	430 / 510	220	380 - 480	14	8
		18															
		20															
		16	16	8 x Q24	505	430 / 510	220	380 - 480	14	8							
		18															
		16	18								8 x Q24	505	430 / 510	220	380 - 480	14	8
		18															
		20															

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)	
116	125	11.1 - 11.6 11.4 - 11.9 11.6 - 12.1	12.3 - 12.8 12.6 - 13.1 12.8 - 13.3	1025 - 1035	600	915 / 995	710 / 780	280	735 / 915	
					650		760 / 830			
				1085 - 1095	600					760 / 830
					650					
					700					
116	125	11.7 - 12.2 12.0 - 12.5 12.2 - 12.7	12.9 - 13.4 13.2 - 13.7 13.4 - 13.9	1055 - 1075	595	920 / 1000	710 / 770	265	765 / 920	
					645		760 / 820			
				1105 - 1125	595					760 / 820
					645					
					695					
116	125	11.0 - 11.5 11.3 - 11.8 11.5 - 12.0	12.2 - 12.7 12.5 - 13.0 12.7 - 13.2	1025 - 1035	600	935 / 1015	750 / 790	280	760 / 930	
					650		800 / 840			
				1085 - 1095	600					800 / 840
					650					
					700					
116	125	8.6 - 9.1 8.9 - 9.4 9.1 - 9.6	13.9 - 14.4 14.2 - 14.7 14.4 - 14.9	1065 - 1070	605	940 / 1020	850 / 1020	280	750	
					655		900 / 1070			
				1115 - 1120	605					900 / 1070
					655					
					705					



KARMA MEDICAL PRODUCTS CO., LTD.
T +886-5-206 6688
Email : globalsales@karma.com.tw

karma
www.KarmaMedical.com