

karma

A Better Wheelchair, A Better Fit



Manual Wheelchair

KM-AT60

Alluring, Stylish and Easy to Carry

With the superior agility and maneuverability design, the KM-AT60 is especially featured with ultra-light weight and extra strength.

This model conquers the difficulties of non-foldable frames of active chair and makes it easy to carry and transport.



KM-AT60

Customized design for seat width, seat height, backrest height, backrest angle, footrest height, camber and adjustable Center of Gravity. Various options for rear wheels and other accessories!



FRAME COLOR



Easy to Carry

Rigid frame with foldable back & quick release rear wheel design facilitates the daily transportation and storage.



The best support

Tension adjustable backrest prevents Hammock Effect and gives users optimal back support.



Multiple choice

Various options for rear wheels and other accessories.

Specifications

WIDTH	F/R wheels												Max. User Weight
16"	5" 2 4Q	38-46	36-51	62	62	N/A	Min.88	46-58	31-48	77	10.5kg	100kg	

KARMA reserves the right to modify information herein without further notice.

±1cm ±0.5kg