

Blazer

One of Karma's longest-running and best-selling product series, Blazer is one of the industry's most compact, agile, and reliable multipurpose power wheelchairs. Whether you're a business professional or an urban explorer, you can choose the seating system, size, and positioning functions you need to make Blazer the perfect everyday wheelchair for you.





Easy Transfer

The flip-up armrests and swing-away footrests make it easy to get on and off the wheelchair and move closer to tables.

Excellent Indoor Maneuverability

One of the most compact power wheelchairs on the market, Blazer handles brilliantly through doorways, into lifts, and around corners.



Outdoor Prowess

Blazer's excellent performance extends outdoors; independent rear wheel suspension and oversized casters ensure a smooth ride, while an optional LED lighting system adds safety for on-road use.

Blazer Sling



Adjustable

Seat depth, back angle, headrest height, headrest angle, armrest height, armpad position, and footrest depth are all adjustable, according to user's needs.



Adaptable

The sling seat is compatible with third-party backrests and accessories.



Compact

Using basic tools, you can quickly and easily fold the sling backrest forward for transport or storage.



Positioning

- Optional Power Tilt-in-Space up to 25°
- Optional Manual Recline 114
- Optional Power Recline 122
- Standard Seat to Back Angle Adjustment 92°/97°/102°



BLAZER

Options and Accessories



Handlebar



Swing-away controller bracket



Ship-shape armpad



Tray table



Elevating Legrests



Push Brake



Push handle for Sling

SIDE FRAME COLORS

- Blue
- Orange
- White

Blazer Captain

Practical

Reflective stripes add safety with style, while the pushbar allows attendants to comfortably maneuver Blazer in freewheel mode.

Adjustable

Armrest height, armpad position, and footrest depth are all adjustable, according to user's needs.

Multiple Sizes

Choose from 450 mm and 475 mm seat widths.



Positioning

Standard Manual Recline up to 134°



Maximize Comfort

Choose either a contoured or flat captain seat with manual recline and width-adjustable armrests.

Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Effective seat width (mm)	Effective seat depth (mm)	Backrest height w/Cushion (mm)	Seat height w/o Cushion (mm)	Seat height w/Cushion (mm)	Power seat lift (mm)	Power Tilt, Seat Plane (X°)	Backrest recline, seat to back angle (X°)	Legrest angle, leg to seat (X°)	F/M/R Wheel (in)	Overall length with legrest (mm)
Blazer Sling	430/480	375-475 425-475	420/470	375-475 425-475	460/540	445	540	-	8-28	92/97/102	105	9/-/14	1090-1160
Blazer Captain	450/475	440/490	470/500	420-500/ 490-570	535	-	565/553	-	-	0-134/ 0-131	100/105	9/-/14	1060-1130

Overall width (mm)	Max. User weight (kg)	Drive range (km)	Max. speed (km)	Obstacle climbing (mm)	Min. turning radius w/o user (mm)	Ability to climb rated slope (X°)	Battery (Ah)	Electrical system
580/600	136	27	10	50	825	12	50*2	VR2/Linx
600/625	136	27	10	50	825	12	50*2	VR2/Linx

- The actual product specifications may differ from the data in the table by ±10 mm ±0.5 kg.
 - The actual product specifications may vary according to different configurations.
 - Karma reserves the right to modify information herein without further notice.

