

# karma

A Better Wheelchair, A Better Fit



Manual Wheelchair

## Champion 100

### *Compact, Comfortable, and Lightweight*

Champion 100 is 30% lighter than a general steel wheelchair. The seamless elliptical-frame design reduces the amount of the welding points, improving the durability and the look of the chair. It is economical, yet it provides multiple ergonomic designs that you will not see on other steel wheelchairs.



# Champion 100

Modern, stylish and unique elliptical design with streamline shape, reliable with low maintenance costs and best choice for your steel wheelchair.

## Ergonomic Armrest Height

The general distance from elbow to hip is around 20-25 cm. However, most of the steel wheelchairs provide 30-35 cm armrest height, which is uncomfortable and may cause shoulder and neck pains. Champion 100 provides a more reasonable armrest height that will make your life easier.



## Seamless Elliptical Frame

The seamless elliptical frame design reduces the amount of the welding points, improving the durability and the look of the chair.

## Optional Attendant Brake



## Ergonomic Push Handle

The handle is designed according to the shape of hands. It is comfortable to hold and easy to push. The attendant brake can be added on upon request.



## Optional Foldable Backrest

The Champion 105 is equipped with foldable backrest, which reduces the overall height and makes it easier to store and to carry around.



## Anti-Slippery Footplate

An anti-slippery footplate provides you a stable and comfortable base for your foot.

## SPECIFICATION

## FRAME COLOR



Model	Width	F/R wheels	Seat Width	Seat Depth	Seat Height	Backrest Height	Backrest Width	Backrest Depth	Backrest Angle	Backrest Tilt	Backrest Lock	Max User Weight
CHM 100	16"/18"	7"/22"	41/46	42	59/62	22	22	94	50	38	86	14 kg / 14.3 kg
CHM 105												14.6 kg / 14.9 kg

±1cm ±0.5kg -Karma reserves the right to modify information herein without further notice. 1111501000127 Rev.00 Mar. 2018



KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

F +886-5-206 7788

Email : globalsales@karma.com.tw

**karma**

www.KarmaMedical.com