

eFlexx

Compact Power Mobility

Adjustable, compact, lightweight, and stylish, eFlexx is the perfect general-use power wheelchair to extend your journey and take you anywhere!



Self-Propel

20" rear wheels with freewheel hubs allow users to easily switch between power drive and self-propelled modes.



Transport

14" rear wheels offer the lightest weight and smallest storage size.



Compact for Travel

In addition to having a foldable frame, eFlexx's battery pack and footrests are removable, making it compact enough to fit into a car.

Quick-Release Battery Pack

The 26 Ah lithium battery provides up to 20 km* of driving range.

The 11.5 Ah Lithium complies with IATA UN38.3 regulations. It allows users to travel conveniently without burden.



Benefits of Lithium Batteries

Lightweight, environmentally-friendly, fast-charging, and long-lasting.

Flexible and Adjustable Design

Configurable seat width and depth, swing-away footrests, backrest angle, and armrest height are all adjustable for a personalized fit.



Swing-away footrests

Enable forward transfers into and out of the wheelchair.



Flip-up armrests

Allow users to laterally transfer into and out of the wheelchair.



Height adjustable armrests

Tool-free height adjustment from 210 – 290 mm to suit user needs.

Foldable push handles

Allow attendants to push the wheelchair and easily fold away when not in use.



Swing-away controller bracket

Approach tables with ease, without damaging your controller or injuring your hand.

User-friendly design

The intuitive, high visibility adjustment levers allow users to easily adjust the armrests and footrests.

OPTIONS

- Wheel Locks
- Rigidifying push bar
- Vertical footrests
- Elevating legrests
- Swing-away controller bracket
- Foldable headrest

ACCESSORIES

- Handle bar
- Cane holder
- IV pole

Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Effective seat width (mm)	Effective seat depth (mm)	Backrest height w/Cushion (mm)	Seat height w/Cushion (mm)	Backrest recline, seat to back angle (X°)	Legrest angle, leg to seat (X°)	F/M/R Wheel (in)	Overall length with legrest (mm)	Overall width (mm)	Max. User weight (kg)
eFlexx-14	360/410/460	415/440/465	405/455/500	560/585/610	410	485	90-110	99	8/-/14	1015-1065	560/615/665	100
eFlexx-20	360/410/460	415/440/465	390/455/510	545/585/610	410	485	90-110	99	8/-/20	1125-1175	590/645/690	100

Drive range (km)	Max. speed (km)	Obstacle climbing (mm)	Min. turning radius w/o user (mm)	Ability to climb rated slope (X°)	Battery (Ah)	Electrical system
10 (11.5AH) 20 (26.1AH)	8	35	795/825/905	10	11.5*1 26.1*1	Linx
10 (11.5AH) 20 (26.1AH)	8	35	795/840/865	10	11.5*1 26.1*1	Linx



Scan to learn more about eFlexx

- The actual product specifications may differ from the data in the table by ± 10 mm ± 0.5 kg.
- The actual product specifications may vary according to different configurations.
- Karma reserves the right to modify information herein without further notice.