



The swingaway footrests and self-propelled rear wheels are detachable, reducing the weight to just 6 kg, making storing and transporting the Ergo Lite 2 a breeze!

Who is **Ergo Lite** for?



Everyday Users

For when a wheelchair user can walk but not far, Ergo Lite's S-Ergo Seating System provides users with all-day comfort and support .



Caregivers

For the carer who cannot carry a lot, Ergo Lite's light weight allows caregivers to preserve their stamina throughout the day.



Travelers

Ergo Lite is the perfect travel companion. Just fold the lightweight wheelchair and put it in the trunk of your car!

NEW

Ergo Lite 2 Tall

The Same Lightweight Freedom, Now with a Better Fit

The Ergo Lite 2 Tall takes our all-time best-selling model to the next level, offering expanded size options, a higher weight capacity, and improved seat-to-floor height. It's the perfect blend of comfort, practicality, and portability for a wider range of users.



Push & Pull Brake

Karma's Push & Pull brake system allows attendants to safely park the wheelchair without bending down to reach the wheel locks.



Crash Tested

The Ergo Lite 2 Tall has been crash tested in accordance with ISO 7176-19, making it suitable for use as a seat in a motor vehicle.




Detachable, Swingaway Footrests

Allows for easy transfer, on and off the wheelchair and detaches with no tools required.



Tool-free quick-release rear wheels. (optional)



With reflective accents to maximize safety on the road.

Easy to Transport and Store

Convenient and suitable for elderly caregivers to take care of elderly users. Easy to carry out and fold for storage.



Scan to learn more about Ergo Lite 2 Tall

Ergo Lite 2

The world's best-selling lightweight travel wheelchair; weighing only 9.2 kg and equipped with detachable, swingaway footrests, the Ergo Lite 2's practicality and portability are best-in-class.



Desk-Length Armrests

Allow the user to get closer to tables.

Push & Pull Brake

Karma's Push & Pull brake system allows attendants to safely park the wheelchair without bending down to reach the wheel locks.



Ergo Handrim

Ergonomic handrims offer improved grip and comfort over traditional designs



Detachable, Swingaway Footrests

Allows for easy transfer, on and off the wheelchair and detaches with no tools required.

Transport

Optional 14" rear wheels for attendant propulsion



Scan to learn more about Ergo Lite 2

Ergo Lite

Karma's lightest wheelchair, weighing a mere 8.6 kg and featuring transport wheels and fixed footrests, Ergo Lite is perfect for those seeking a durable travel wheelchair, and because its crash test approved, you can go further than ever!



Should have done this sooner..

" I am a caregiver rather than a patient; this chair has very likely saved my back. The chair that insurance provided was a whopping 42 lbs; I am 60, transporting my post-stroke husband as we try to keep a reasonably active schedule. He's gaining strength, but balance and vision issues make a chair a must. I was managing my discomfort, but knew as it got worse that a change was necessary. "

Lovisa P / United States



Anne Crofts / U.K.



Arma Hong / Thailand



A really great chair !!!

" I have had 5 wheelchairs that have ranged in price from about \$200 to well over \$1000, and this is by far the nicest. It is built well, has large and very comfortable armrests, and is just great looking!! The very best features are that it is about 3" narrower than my previous chairs (less damage to walls and door frames) and it is 14 lbs lighter than my last chair. I love it !!! "

quadman / United States



Scan to learn more about Ergo Lite

Series	Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)
Ergo Lite	Ergo Lite 2 Tall	410	435	7 / 14	510	435	210	410 - 460	15	6
		460								
		510								
		410		7 / 22						
		460								
		510								
	Ergo Lite 2	400	420	6 / 14	490	420	210	360 - 410	11	6
		450		6 / 20						
		400								
		450								
Ergo Lite	400	420	6 / 14	490	420	210	395 - 420	11	6	
	450									

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)
100	115	6.9 (14Q) 8.7 (14F)	9.9 (Q) 9.8 (F)	990 - 1005	565	950	720 (Q) 647 (F)	230 (Q) 240 (F)	690 (Q) 745 (F)
		7.1 (14Q) 8.8 (14F)	10.1 (Q) 9.9 (F)	1015 - 1030	615				725 (Q) 745 (F)
		7.4 (14Q) 9.2 (14F)	10.4 (Q) 10.3 (F)		665			260 (Q) 270 (F)	760 (Q) 785 (F)
		7.1 (22Q) 10.1 (22F)	11.2 (Q) 11.1 (F)	990 - 1005	575		720 (Q) 755 (F)	260 (Q) 280 (F)	690 (Q) 745 (F)
		7.2 (22Q) 10.2 (22F)	11.4 (Q) 11.3 (F)	1015 - 1030	625				725 (Q) 745 (F)
		7.7 (22Q) 10.7 (22F)	11.9 (Q) 11.8 (F)		700			760 (Q) 785 (F)	
100	100	6.7 (14Q) 8.1 (14F)	9.2	962 - 970	570	900	635	280	700
		6.9 (14Q) 8.3 (14F)	9.3		615			295	710
		6.7 (20Q) 9.7 (20F)	10.7	960 - 975	595		720	280	700
		6.9 (20Q) 9.8 (20F)	10.9		645			295	710
100	100	8.6	8.6	962 - 970	560	905	852	245	710
		8.8	8.8	962 - 970	610				

Rev.04 Jun . 2025

KARMA MEDICAL PRODUCTS CO., LTD.
T +886-5-206 6688
Email : globalsales@karma.com.tw

karma
www.KarmaMedical.com