

# Ergo Live

## Live Active in Style and Comfort

Ergo Live is highly configurable to provide the ideal balance between comfort and performance. The lightweight frame is foldable and features a minimalist design with a splash of color, allowing you to go anywhere in style.

## Ergonomic Seating



### 1 S-Ergo Seating

Evenly distribute users' weight from the hip bone to the lower thigh for comfort and support.



### 2 AEGIS Microbe Shield

Washable and breathable upholstery that inhibits bacterial growth to eliminate odor and protect your health.



### 3 Tension Adjustable Backrest

Series of individually-adjustable Velcro straps create a user-specific fit to maximize comfort and support.



### 4 Adjustable Footplate Angles

The one-piece footplate is height- and angle-adjustable to offer the support you need.



## Configuration




### 1 Adjustable Center of Gravity Active / Neutral / Stable



### 2 Adjustable Backrest Angle 3-stage adjustment -5° / 0° / 5°



### 3 Adjustable Backrest Height 360 - 455 mm



### 4 Adjustable Armrest Height 225 / 250 mm



### 5 Configurable Seat Height 445 - 510 mm (with different caster fork)



### 3 Multiple Seat Dimensions





## Optional Accessories




5" Aluminum Rim Casters



24" Lively Rear Wheel  
(24-spoke)



Push Handles



Foldable Push Handles



Dual Function  
Skirt/Mud Guard



Flip-up Armrest



2D Flip-up One Piece  
Footplate

## FRAME COLOR

- White
- Silver Black
- Lemon Yellow



Scan to learn more  
about Ergo Live

Series	Product Name	Nominal seat width (mm/in)	Nominal seat depth (mm/in)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)
Active & Specialty	Ergo Live	340 / 14	390 / 15	5 x 24	465	360 - 455	225 / 250	375 - 435	7 / 12 / 17	5
		340 / 14	415 / 16							
		390 / 16	415 / 16							
		390 / 16	465 / 18							
		440 / 18	415 / 16							
		440 / 18	465 / 18							

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)
97	125	9.0	13.1	895 - 910	580	810 - 885	855	290	550
		9.1	13.2	920 - 935			870		
		9.1	13.3	920 - 935	630		870	310	585
		9.2	13.5	970 - 985			900		
		9.2	13.5	920 - 935	680		870	330	620
		9.3	13.7	970 - 985			900		

Rev.04 Jun . 2025