

Flexx Junior



Adjustable Seat Dimensions

Adjustability plays a key role in pediatric wheelchairs; the seat width of the Flexx Junior can be set between 255 - 355 mm, while the seat depth can be set from 285 - 365 mm.



The standard tension adjustable backrest can be configured according to each user's needs. The seat to back angle is adjustable from -4 to 16°.



90° footrests and adjustable footplates can be set to different positions and angles.



Easy Transfers

Flip-up armrests and detachable, swingaway footrests improve convenience and easier access to the chair.

Other Adjustments Inherited from Flexx

The seat height is adjustable between 394 - 457 mm, depending on the caster and rear wheel size and position, and the backrest angle can be set between -4° and 16°.



Transport

Attendant-propelled with 14" rear wheels

OPTIONS & ACCESSORIES

- Lap Belt
- Pelvic Belt
- Transparent Tray Table
- I.V. Pole
- Oxygen Tank Holder
- Cane Holder
- Wood Tray Table (330 and 355 mm seat widths only)
- Push Handle Extension
- Height Adjustable Backrest
- Elevating Detachable Legrest
- Height Adjustable and Detachable Armrest



Multi-functional to Fit Children's Needs

The Flexx Junior is a multi-adjustable pediatric wheelchair that adapts with the child as they change and grow. The adjustable seat width, depth, and height are just a few of the many built-in features that this chair offers. In addition, the Flexx family's accessories are compatible with the Flexx Junior.

SPECIFICATION

Nominal seat width (mm)*	255 - 355
Nominal seat depth (mm)**	285 - 365
Front caster / Rear wheel size (in)	5/6 x 14 , 5/6 x 20 , 5/6 x 22
Front seat height (mm)	410 , 400 , 415
Backrest height (mm)	380
Armrest to seat distance (mm)	155 - 235 / 180 - 235 / 220 - 235
Footrest to seat distance (mm)	150 - 250
Max. user weight (kg)	75
Mass of the heaviest part (kg)	9.7 - 9.8 / 9.4 - 9.5 / 9.4 - 9.5
Total mass (kg)	14.3 - 14.4 / 14.9 - 15.0 / 15.2 - 15.3
Overall L x W x H (mm)	Min. 820 x 450 x 800 / Max. 890 x 590 x 875
Overall folded L x W x H (mm)	Min. 810 x 325 x 725 / Max. 870 x 355 x 725

* Nominal seat width is measured the width of seat cushion at 120mm in front of the backrest. This value is mainly used in the product ordering process.


** Nominal seat depth is measured from the front end of the seat cushion to the back cushion. This value is mainly used in the product ordering process.

- The actual product specifications may differ from the data in the table by ± 1 cm ± 0.5 kg.

- The actual product specifications may vary according to different configurations.

- Karma reserves the right to modify information herein without further notice.

FRAME COLORS



Scan to learn more about Flexx Junior

Rev.04 Jun . 2025

KARMA MEDICAL PRODUCTS CO., LTD.
T +886-5-206 6688
Email : globalsales@karma.com.tw

karma
www.KarmaMedical.com