

# KM-1500



## Standard Version and Easy to Store

KM-1500 is an ideal choice for short term users, healthcare centers, and hospitals. With a lightweight aluminum frame and foldable components, the KM-1500 caters to basic mobility needs.



### Foldable Backrest

Reduces storage space and makes it easy to lift.



### Swivel Footplates

Provide more space for leg movement and allows for easy transfer, on and off the wheelchair.

### Attendant Brakes

Dual brake system brings more convenience for caregivers and enhances the safety of the user.



## SPECIFICATION

Nominal seat width (mm)*	410 / 460
Nominal seat depth (mm)**	415
Front caster / Rear wheel size (in)	8 x 24
Front seat height (mm)	485
Backrest height (mm)	380
Armrest to seat distance (mm)	225
Footrest to seat distance (mm)	295 - 415
Max. user weight (kg)	100
Mass of the heaviest part (kg)	13.8 / 13.9
Total mass (kg)	14.1 / 14.2
Overall L x W x H (mm)	Min. 1010 x 720 x 865 / Max. 1060 x 770 x 865
Overall folded L x W x H (mm)	915 x 350 x 770

\* Nominal seat width is measured the width of seat cushion at 120mm in front of the backrest. This value is mainly used in the product ordering process.

\*\* Nominal seat depth is measured from the front end of the seat cushion to the back cushion. This value is mainly used in the product ordering process.

- The actual product specifications may differ from the data in the table by  $\pm 1$  cm  $\pm 0.5$  kg.

- The actual product specifications may vary according to different configurations.

- Karma reserves the right to modify information herein without further notice.



Scan to learn more about KM-1500

Rev.04 Jun . 2025

KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

Email : globalsales@karma.com.tw

**karma**  
www.KarmaMedical.com