

# KM-2500

## Ultra Narrow and Easy to Transport

Weighing 10.9 kg, the KM-2500 is an ideal choice for people with a petite frame or for people living in tight spaces with narrow doorways. It features a foldable backrest and swing away footplates. Made of aluminum alloy, it is lightweight and compact, making it easy to store and transport.



### Outrigged Armrest

The outrigged armrests minimize the chair's overall width while still maintaining an ergonomic sitting space.



### Self-Propel Option

Optional 20" rear wheel with Ergo Handrim for self-propelling.



### Foldable Backrest

Compact and foldable backrest reduces storage space and makes it easy for traveling.

## Clean, Comfortable, and Fresh

- Anti-microbial upholstery inhibits bacteria, fungus, and eliminates odor.
- Removable & washable
- Breathable



## SPECIFICATION

Nominal seat width (mm)*	420 / 470
Nominal seat depth (mm)**	400
Front caster / Rear wheel size (in)	6 x 14
Front seat height (mm)	450
Backrest height (mm)	420
Armrest to seat distance (mm)	225
Footrest to seat distance (mm)	310 - 410
Max. user weight (kg)	100
Mass of the heaviest part (kg)	10.7 / 10.8
Total mass (kg)	11.2 / 11.3
Overall L x W x H (mm)	Min. 885 x 550 x 855 / Max. 915 x 600 x 855
Overall folded L x W x H (mm)	740 x 250 x 680



Scan to learn more  
about KM-2500

\* Nominal seat width is measured the width of seat cushion at 120mm in front of the backrest.  
This value is mainly used in the product ordering process.

\*\* Nominal seat depth is measured from the front end of the seat cushion to the back cushion.  
This value is mainly used in the product ordering process.

- The actual product specifications may differ from the data in the table by  $\pm 1\text{ cm} \pm 0.5\text{ kg}$ .

- The actual product specifications may vary according to different configurations.

- Karma reserves the right to modify information herein without further notice.

Rev.04 Jun . 2025

KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

Email : globalsales@karma.com.tw

**karma**  
www.KarmaMedical.com