

KM-5000



Standard Reclining Wheelchair

The KM-5000 is a standard version of the recline series, designed for long-term users; fully detachable armrests allow for lateral transfers.



Detachable Armrest

Fully detachable armrests offer greater convenience when transferring in and out of the wheelchair.



Elevating Legrests

The leg rest angle can be adjusted for better leg support when reclining.



Clean, Comfortable, and Fresh

- Anti-microbial upholstery inhibits bacteria, fungus, and eliminates odor.
- Removable & washable
- Breathable

SPECIFICATION

Nominal seat width (mm)*	425 , 460
Nominal seat depth (mm)**	420 , 470
Front caster / Rear wheel size (in)	8 x 14 / 8 x 24
Front seat height (mm)	500
Backrest height (mm)	450
Armrest to seat distance (mm)	235 / 230
Footrest to seat distance (mm)	450 - 550 / 375 - 425
Max. user weight (kg)	115
Mass of the heaviest part (kg)	13.2 / 13.5 / 15.7 / 16.0
Total mass (kg)	18.7 / 19.0 / 21.3 / 21.6
Overall L x W x H (mm)	Min. 1215 x 660 x 1295 / Max. 1290 x 660 x 1295
Overall folded L x W x H (mm)	Min. 745 x 330 x 935 / Max. 960 x 400 x 930

* Nominal seat width is measured the width of seat cushion at 120mm in front of the backrest. This value is mainly used in the product ordering process.

** Nominal seat depth is measured from the front end of the seat cushion to the back cushion. This value is mainly used in the product ordering process.

- The actual product specifications may differ from the data in the table by ± 1 cm ± 0.5 kg.
- The actual product specifications may vary according to different configurations.
- Karma reserves the right to modify information herein without further notice.



Scan to learn more about KM-5000

Rev.04 Jun . 2025

KARMA MEDICAL PRODUCTS CO., LTD.
T +886-5-206 6688
Email : globalsales@karma.com.tw

karma
www.KarmaMedical.com