



KARMA MEDICAL PRODUCTS CO., LTD.

NO.2363 Sec,2 University Rd.,Min-Hsiung Shiang,Chia-Yi 621,Taiwan

T +886-5-206 6688 F +886-5-206 7788

Email: globalsales@karma.com.tw
www.KarmaMedical.com



Karma reserves the right to modify information herein without further notice 104-40600 Rev. 01 May. 2015

karma A Better Wheelchair, A Better Fit



Power Wheelchair

KP-25.2

*Discover the path and make it yours
Easy access with longer driving range*

KP-25.2 is armrest combined with ergonomically designed arm-pad enables flipping and detachments, offering optimal comfort according to your needs.

It is ultra-lightweight, foldable, and has easy battery removal, making transport even more convenient.





◀ **Explore the urban jungle in a compact mid-size wheelchair**

The overall width of the wheelchair is 59 cm.



▲ **" Let's go! "**

The heaviest part weighs less than 33kg, for easy transport and storage.



																±1cm ±0.5kg	
Seat Width	Seat Height	Overall Length	Overall Width	Overall Height		MOTOR	BATTERY	WHEELS (F/M/R)		Turning Radius	Max Safe Slope		KM / HR	RANGE	MAX	Weight with Batteries	Weight w/o Batteries
42	50cm	110cm	59cm	96cm	39~51cm	200W	50Ah	8"12"	VR2	82cm	10°	5cm	6Km	38Km	115kg 135kg(Pneumatic)	70kg	35kg
47	50cm	110cm	64cm	96cm	39~51cm	200W	50Ah	8"12"	VR2	83cm	10°	5cm	6Km	38Km	115kg 135kg(Pneumatic)	70kg	35kg

• KARMA reserves the right to modify information herein without further notice.

Features

- Aluminum alloy frame prevents rusting
- Easy battey access mechanism
- Foldable backrest
- Multi-functional controller
- Multi-dimension adjustable armrest



Flip-back armrest



Height-adjustable armrest



Easy battery removal



Detachable Footrest



Aluminum alloy split-rims



Easy access free wheel levers



Tool free quick folding for easy transport

Optional accessories

- **Seat Pan**
For aftermarket cushions



- Quick release attendant controller
Allows the joystick to be moved from front to back in seconds.



- Swing away controller bracket



- Adjustable headrest



- Curb climber
enables crossing one single step (7-10cm)



- Foldable drink holder



- Amputee footrest



- Elevating footrest