

KS-200

The most convenient and transportable mobility scooter.



Shopping, Travelling, Living a Happy and Carefree Life

Let's go !



An achievement in small mobility scooter category. With ultra-compact size and newly developed quick-release design, KS-200 creates the best way to travel in an urban setting!

3 Easy steps to disassemble

1



Removable seat

Hold the seat and pull the handle below, lift it up to separate.

2



Detachable battery pack

Rotate the switch counterclockwise, pull up to release the battery pack.

Charge up easily

It's able to recharge the battery pack individually.

3



Easy to Dismantle

Pull up the red handle to separate the front and rear frame.



Lightweight and Compact Design

- The heaviest piece is only 13.3kgs after disassembly.
- Optimal size for easy to go in and out of the elevator.



Overall Width	Overall Length	Overall Height	F/R Wheels	Motor	Battery	Charger	Min. Turning Radius	Max. Safe Slope	Max. Speed	Range	Max. User Weight	Weight	Weight w/o Battery
55cm	104cm	92cm	9" / 9"	270W	22Ahx2	2Ahp	147.5cm	8°	6KM/H	15KM/H	115kg	51.5kg	38.2kg

Karma reserves the right to modify information herein without further notice.



ISO 9001

karma

E-mail:globalsales@karma.com.tw
www.KarmaMedical.com