

Create ultimate experience from the heart. Through the distinguished mechanical design, we make this mid-wheel-drive wheelchair with the best maneuverability and high performance. Furthermore, its full-scale adjustable seating system with power position changes will provide you the most comfortable experience than ever. Perform a wonderful spin with MID Lectus!



LECTUS Seat

Release body pressure with the Lectus positioning seating system, its multiple adjustments and smooth texture support a pleasurable sitting.



MID Chassis

With silent chassis and 6 wheel linkage suspension system, MID Lectus brings you the utmost comfort, stable and smooth maneuverability for both indoor and outdoor use.

6 WHEEL SUSPENSION SYSTEM



Designed with linked suspension on the front and rear caster, it can evenly distribute the driving impact and enable to overcome curbs up to 7.5 cm.



45° Power Tilt-in-space

Release the pressure at buttocks and stabilize the body position while driving downhill.



30cm Power Elevating System

Easy to reach higher spaces.



180° Power Recline + 180° Power Elevating Legrest

Total flat position allows the user to lie down straight.



Advanced Driving Range

Larger battery box that can carry 80Ah battery, with embedded battery design provides tool-free assembly for convenient maintenance. Go farther from now on!

OPTIONS & ACCESSORIES

- Transfer Handles
- Push Handles
- Rear Bumper Roller
- Parallel Swing-Away-Down Controller Bracket
- Table Width Flip Around Control
- 3D Adjustable Headrest
- 3D Lateral Support
- 3D Hip Support
- 3D Calf Support
- Positioning Belt
- Positioning Belt Retractable
- LED Lighting System Include Indicators

Special controls & individual adoptions are available on request

SPECIFICATION

Seat Width (cm)	42 - 52
Front/Middle/Rear Wheel (inch)	8 / 14 / 6
	Tilt 0°-45°
Power Functions	Bio-mechanical Recline (-5°)-90° Bio-mechanical Legrest: (-5°)-90° Lift: 0-30cm
Controller	Rnet
Motor	350W
Battery	50Ah*2 / 80Ah*2
Min. Turning Radius (cm)	62
Min. Ground Clearance (cm)	9
Max. Safe Slope	10°
Max. Obstacle Height (cm)	7.5
Max. Speed (km/h)	6,10,12
Max. Driving Range (km)	24/36
Max. User Weight (kg)	140
w/ Batteries / w/o Batteries (kg)	190(80Ah)/ 133
Overall Length x Width x Height (cm)	110*61.5*104

±1cm ±0.5kg  
-Karma reserves the right to modify information herein without further notice.  
-For more details, please refer to the Specification Pages No. 36-37.  
-Battery is not part of the product, might change as per local country specification and availability.