

# S-Ergo 125

## Comfortable Seating and Easy Transfer

The classic model of the S-Ergo line, the S-Ergo 125 features flip-up armrests and detachable, swing-away footrests, providing users with enhanced comfort and ease during daily activities and transfers.



### AEGIS Microbe Shield

Named for Zeus' mythical shield, AEGIS Microbe Shield upholstery treatment inhibits bacterial growth to eliminate odor and protect your health.



### Ergo Armpads

The ergonomic armpads are curved inward according to the way users lean their forearms and elbow on the pad.



The armrests flip up for easier transfer in and out of the wheelchair.



### Ergo Handrim

The ergonomic handrims improve propelling comfort and efficiency.



### Detachable, Swingaway Footrests

Enable easy transfer, on and off the wheelchair and detaches with no tools required.



### Compact

Lightweight & compact, easily folded, convenient for storage and transportation.



## S-Ergo 115

Fixed armrests and detachable, swingaway footrests make this an ideal option for forward transfers.



## S-Ergo 105

The most affordable S-Ergo wheelchair, featuring fixed armrests and fixed footrests.



Series	Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)				
S-Ergo	S-Ergo 125	400	430	7 / 14	500	430	230	335 - 415	12	6				
		450												
		500												
		400												
		450												
		500												
		440	470	7 / 14	515	425		360 - 435	14	8				
		500												
		440		7 / 24					515	425	360 - 435	13	9	
		500												
	500													
	S-Ergo 115	400	430	7 / 14	500	430	220	355 - 430	12	6				
		450												
		500												
		400	420	7 / 24	495						355 - 430	12	6	
		450												
		500												
		440	470	7 / 14	515	425		230	360 - 435	14	8			
500														
440		7 / 24		515						425	230	360 - 435	14	9
500														
500														

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)		
97	115	9.1 (14Q) 11.3 (14F)	12.9	1010	575	905	735	310	705		
		9.3 (14Q) 11.5 (14F)	13.1		625						
		9.5 (14Q) 11.7 (14F)	13.3		675						
		95	115	9.1 (24Q) 13.0 (24F)	14.6	1000	580	910	805	320	710
				9.4 (24Q) 13.2 (24F)	14.8		630				730
				9.6 (24Q) 13.4 (24F)	15.0	680	750				
10.2 (14Q) 12.7 (14F)	14.5			1090	630	975	845				350
10.3 (14Q) 13.0 (14F)	14.8	680	780								
9.7 (24Q) 13.5 (24F)	15.3	655	775								
9.9 (24Q) 13.7 (24F)	15.5	705	800								
95	115	8.1 (14Q) 10.3 (14F)	11.9	1010	575	905	735	310	705		
		8.3 (14Q) 10.5 (14F)	12.1		625						
		8.5 (14Q) 10.7 (14F)	12.3		675						
		95	115	8.1 (24Q) 12.0 (24F)	13.6	1000	580	910	805	320	710
				8.4 (24Q) 12.2 (24F)	13.8		630				730
				8.6 (24Q) 12.4 (24F)	14.0	680	750				
95	115			9.3 (14Q) 11.7 (14F)	13.5	1090	615				975
		9.5 (14Q) 11.9 (14F)	13.7	665	780						
		9.1 (24Q) 13.0 (24F)	14.8	640	760						
		9.3 (24Q) 13.2 (24F)	15.0	690	800						

Series	Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)
S-Ergo	S-Ergo 105	400	430	7 / 14	500	430	230	360 - 435	13	7
		450								
		400		7 / 24						
		450								

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)
98	115	11.5	12.1	1010	585	900	910	325	735
		11.7	12.3		635				
		9.4 (24Q) 13.2 (24F)	13.8	1005	590	905	980	335	740
		9.6 (24Q) 13.4 (24F)	14.0		640				

Rev.04 Jun . 2025



Scan to learn more about S-Ergo 125



Scan to learn more about S-Ergo 115



Scan to learn more about S-Ergo 105

KARMA MEDICAL PRODUCTS CO., LTD.  
T +886-5-206 6688  
Email : globalsales@karma.com.tw

**karma**  
www.KarmaMedical.com