

# karma

A Better Wheelchair, A Better Fit



Manual Wheelchair

## S-Ergo 205

### *Compact, Comfortable, and Lightweight*

S-Ergo 205 is designed for those who seek easy transport in everyday use. Worldwide patent S-Ergo seating system provides excellent pressure relief at the hipbone area. The outriggered armrest design creates a comfortable sitting space with the most compact total width possible.



# S-Ergo 205

Two rear wheel sizes upon request.  
14" for easy transit  
20" for self-propel



ERG 205-14

ERG 205-20



### S-Ergo seating system, good pressure relief

Worldwide patented S-Ergo shaped seating system will allow pressure to be redistributed and help the user to stay in position instead of sliding down the seat, and is pressure mapped for the best possible shape.



### Compact and foldable, easy for care-givers

The aluminum alloy reduces the weight significantly. Foldable backrest makes the wheelchair more convenient to transport and store.



### Outrigged armrest, more space, and more comfort

Outrigged design minimizes the chair's overall width while maintaining an ergonomic and comfortable sitting space.



### Foldable cup holder (Optional)

Meet your daily need without increasing the storage space.




### Aegis microbe shield

The upholstery is treated with Aegis microbe shield, which acts as an anti-bacterial shield for improved hygiene.



### Outward swivel footplates

Easier for users to transfer from the wheelchair to other place.



### Brilliant armrest curve design

Enabling the user to get closer to the environment and providing additional space for the thigh.

## SPECIFICATION

	F/R wheels	Seat width	Seat depth	Seat height	Backrest height	Wheelchair width	Wheelchair depth	Wheelchair height	Seat height	Seat depth	Seat weight	Max. User Weight
16"	6"x20"	41	42	57	31	22	93	47	39	85	10.8 kg	100 kg
18"	6"x20"	46	42	62	31	22	93	47	39	85	11 kg	100 kg
16"	6"x14"	41	42	57	31	22	89	47	39	85	9.8 kg	100 kg
18"	6"x14"	46	42	62	31	22	89	47	39	85	10 kg	100 kg