

# karma

A Better Wheelchair, A Better Fit



Manual Wheelchair

## S-Ergo 305

### *Comfort, Adjustability, and Ruggedness.*

Sturdy yet lightweight, a premium choice for self-propel wheelchair. S-Ergo seating system helps to redistribute pressure and promote a better seating position.

Adjustable flip-back armrests, swing-away footrests, and quick-release rear wheels are just a few amazing designs that will make your experience more comfortable.



# S-Ergo 305

## Patented S-Ergo seating system. Great pressure relief.

Worldwide patented S-Ergo shaped seating system help the user to stay in position instead of slumping forward. It redistributes the sitting pressure from the hipbone area to the entire thigh, reducing the risk of developing pressure sores.



### Height adjustable and flip-back armrests

A suitable armrest height protects user from shoulder pain. The flip-back armrest opens up space for user when transferring.



### Swing in/out footrest

Easy to swing inward, outward, and to detach, offering smoother transferring in confined environment.



*"The S curved seat is odd at first, but so very very comfortable to her. She no longer complains of back pain from sitting in the chair for long periods."*

*-kindle customer, amazon.com*



### Brand new mudguard

Prevents user from getting dirty when propelling on the road.



### Multiple adjustable mechanism

16" & 18" seat widths available. Armrest height, seat height, and footrest height are all adjustable to fit individuals.

### ● Quick release rear wheels

Reduce the total weight and size of the chair. Easy to carry and store. Ergo Handrim is easy to grip and propel.

### ● AEGIS microbe shield

The upholstery is treated with Aegis microbe shield, which acts as an anti-bacterial shield for improved hygiene.

## SPECIFICATION

	F/R wheels											Max. User Weight
16"	6"x24"Q	41	43	57	33	22-27	108	51.5	40	97	14.2 kg	115 kg
16"	7"x24"Q	41	43	57	33	22-27	108	49.5	40	97	14.2 kg	115 kg
18"	6"x24"Q	46	43	63	33	22-27	108	51.5	40	97	14.6 kg	115 kg
18"	7"x24"Q	46	43	63	33	22-27	108	49.5	40	97	14.6 kg	115 kg

±1cm ±0.5kg -Karma reserves the right to modify information herein without further notice. 111501000132 Rev.00 Dec. 2017



KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

F +886-5-206 7788

Email : globalsales@karma.com.tw

# karma

www.KarmaMedical.com