

SOMA 103

Lightweight aluminum wheelchair featuring fixed desk-length armrests and fixed footrests with flip-up footplates.



Detachable Seating System
Innovative Detachable Seating
Intuitive design makes it easy to detach and clean the upholstery.



Foldable Backrest
For easy transport and storage.



Anti-Slip Footplates
Allow users to effortlessly maintain proper foot positioning throughout the day.



Desk-length Armrest Design
Allows the wheelchair user to roll closer to a desk or table.



SPECIFICATION

Nominal seat width (mm)*	420
Nominal seat depth (mm)**	430
Front caster / Rear wheel size (in)	7 x 16 / 7 x 22
Front seat height (mm)	500
Backrest height (mm)	400
Armrest to seat distance (mm)	235
Footrest to seat distance (mm)	395 - 420
Max. user weight (kg)	100
Mass of the heaviest part (kg)	10.4 / 11.6
Total mass (kg)	10.4 / 11.6
Overall L x W x H (mm)	945 x 585 x 875 / 980 x 615 x 865
Overall folded L x W x H (mm)	945 x 280 x 865 / 940 x 295 x 865

* Nominal seat width is measured the width of seat cushion at 120mm in front of the backrest. This value is mainly used in the product ordering process.

** Nominal seat depth is measured from the front end of the seat cushion to the back cushion. This value is mainly used in the product ordering process.

- The actual product specifications may differ from the data in the table by ± 1 cm ± 0.5 kg.
- The actual product specifications may vary according to different configurations.
- Karma reserves the right to modify information herein without further notice.

16" transport rear wheel option



Scan to learn more about SOMA 103

Rev.04 Jun . 2025

KARMA MEDICAL PRODUCTS CO., LTD.
T +886-5-206 6688
Email : globalsales@karma.com.tw

karma
www.KarmaMedical.com