



Manual Wheelchair



Standing wheelchair

The ability to stand up strongly promotes user independence, increases health benefits and improves life imbalances.

The user will feel comfortable and confident in a near-natural standing position.



BENEFITS

- Activates circulation system and improves cardio-pulmonary function.
- Lowers the risk of developing decubitus sores.
- Improves bowel regularity and function.
- Improves bladder and intestinal motility.
- Reduces contractures and ankylosis.
- Prevents atrophy of the leg muscles.

FEATURES

- 1) Manual wheelchair with power control for standing up.
- 2) Four-point touch ground in standing position providing extra safety.
- 3) Ultra-light weight.
- 4) 80% of the frame is made of ultra light aircraft-grade aluminum.
- 5) 25.8 kg without batteries.
- 6) Adjustable seat depth from 44 cm, 47 cm to 50 cm.
- 7) Seat width: 40 cm, 44 cm or 48 cm.
- 8) Ergonomic flip-back armrest with height-adjustable range from 18 cm to 26 cm.
- 9) Foldable footrest with angle adjustment characteristics.
- 10) Angle & depth adjustable knee pad to fit individual preferences.



± 1 cm ± 0.5 kg

F/R wheels													Max. User Weight
6"x24"	40/44/48	44/47/50	59/63/67	18~26	98.5~104.5	52	42	96/143	2.7Ah x 2	28.6 kg	25.8 kg without batteries		120 kg

KARMA reserves the right to modify information herein without prior notice.

CE   ISO 9001