

SOMA 105

Lightweight aluminum wheelchair featuring fixed full-length armrests and fixed footrests with flip-up footplates



Foldable Backrest
For easy transport and storage.



Ergonomic Push Handle
Maximizes attendant comfort when maneuvering the chair and braking.



Attendant Tipping Lever
Provides additional leverage to easily overcome obstacles.



Modern, stylish, and patented elliptical design.



SPECIFICATION

Nominal seat width (mm)*	400 , 450
Nominal seat depth (mm)**	410 / 440
Front caster / Rear wheel size (in)	7 x 16 / 7 x 22
Front seat height (mm)	500 / 480
Backrest height (mm)	385 / 390
Armrest to seat distance (mm)	220
Footrest to seat distance (mm)	460 - 510 / 395 - 445
Max. user weight (kg)	100
Mass of the heaviest part (kg)	10.3 / 10.6 / 11.2 / 11.5
Total mass (kg)	10.5 / 10.8 / 11.4 / 11.7
Overall L x W x H (mm)	Min. 920 x 595 x 880 / Max. 975 x 665 x 875
Overall folded L x W x H (mm)	840 x 300 x 775 / 955 x 285 x 760

16" transport rear wheel option



Scan to learn more about SOMA 105

* Nominal seat width is measured the width of seat cushion at 120mm in front of the backrest. This value is mainly used in the product ordering process.
 ** Nominal seat depth is measured from the front end of the seat cushion to the back cushion. This value is mainly used in the product ordering process.
 - The actual product specifications may differ from the data in the table by ± 1 cm ± 0.5 kg.
 - The actual product specifications may vary according to different configurations.
 - Karma reserves the right to modify information herein without further notice.

Rev.04 Jun . 2025

KARMA MEDICAL PRODUCTS CO., LTD.
 T +886-5-206 6688
 Email : globalsales@karma.com.tw

karma
 www.KarmaMedical.com