



Power Wheelchair

SP-100

The Easiest Power Wheelchair Solution

SP-100 is a foldable power wheelchair designed for those who look for a simple, lightweight, yet cost-effective solution.

Within few seconds, you can remove the batteries, put the chair into your trunk, and head out for a great day. SP-100 is designed to be light and easy-to-travel without compromising its durability.



SP-100

Best Choice for City Travelers

SP-100 is one of the most compact batteries-removable power wheelchairs on the market. It can even be folded with the batteries loaded, which makes it the perfect choice for the travelers.



Removable battery

22AH*2 Battery Capacity-
Range 15 km.



Elliptical tubes

Lightweight Robust-Structure Frame

The innovative elliptical tubes make the chair light and durable at the same time.



Reinforced Double Cross-Braces Design

SP-100 has completed Two-Drum Test for 200 thousands rounds and 6666 times Drop Test, ensuring the extra durability and user safeness.



Armrests

Easy Transfer

The flip-back armrests provide easy transfer to the chair.



Get Ready to Travel Wherever and Whenever You Want

Linx Controller

5-stage adjustable driving speed ,
with the max. speed of 6km/hr.

Solid casters and wheels — flat-free

















The SP-100 is designed to handle all the city terrain with flat-free 8" caster and 16" rear wheel.



Comfortable and Easy to Clean

Double-layered seat upholstery
is removable and easy to clean.

SPECIFICATION

																
41/46cm	42cm	49cm	99cm	56/61cm	86cm	37-43cm	100 W	22AH*2	8"/16"	Linx	81 cm	6°	4 cm	6	15km	100kg
Weight with Batteries	Weight w/o Batteries															
43kg	29kg															

+1cm ±0.5kg. Karma reserves the right to modify information herein without further notice.

±1cm ±0.5kg -Karma reserves the right to modify information herein without further notice.



KARMA MEDICAL PRODUCTS CO., LTD.

T +886-5-206 6688

F +886-5-206 7788

Email : globalsales@karma.com.tw

karma

www.KarmaMedical.com