

Series	Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)
Ergo Lite	Ergo Lite 2 Tall	410	435	7 / 14	510	435	210	410 - 460	15	6
		460								
		510								
		410		7 / 22						
		460								
		510								
	Ergo Lite 2	400	420	6 / 14	490	420	210	360 - 410	11	6
		450		6 / 20						
		400								
		450								
	Ergo Lite	400	420	6 / 14	490	420	210	395 - 420	11	6
		450								
S-Ergo	Ergo 3	410	435	6 / 14 7 / 14	435 485	435	220	385 - 460	20	6
		460								
		510								
		410		6 / 20 7 / 20	435 485					
		460								
		510								
		410		7 / 24	485					
		460								
		510								
	Ergo 2	410	435	6 / 14 7 / 14	435 485	435	220	385 - 460	20	6
		460								
		510								
		410		6 / 20 7 / 20	435 485					
		460								
		510								
		410		7 / 24	485					
		460								
		510								
	S-Ergo 125	400	430	7 / 14	500	430	230	335 - 415	12	6
		450								
		500								
		400		7 / 24	495					
		450								
		500								
		440	470	7 / 14	515	425		360 - 435	14	8
		500								
		440		7 / 24					13	9
		500								
		500								
	S-Ergo 115	400	430	7 / 14	500	430	220	355 - 430	12	6
		450								
		500								
		400	420	7 / 24						
		450								
		500								

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)
100	115	6.9 (14Q) 8.7 (14F)	9.9 (Q) 9.8 (F)	990 - 1005	565	950	720 (Q) 647 (F)	230 (Q) 240 (F)	690 (Q) 745 (F)
		7.1 (14Q) 8.8 (14F)	10.1 (Q) 9.9 (F)	1015 - 1030	615				725 (Q) 745 (F)
		7.4 (14Q) 9.2 (14F)	10.4 (Q) 10.3 (F)		665			260 (Q) 270 (F)	760 (Q) 785 (F)
		7.1 (22Q) 10.1 (22F)	11.2 (Q) 11.1 (F)	990 - 1005	575		720 (Q) 755 (F)	260 (Q) 280 (F)	690 (Q) 745 (F)
		7.2 (22Q) 10.2 (22F)	11.4 (Q) 11.3 (F)	1015 - 1030	625				725 (Q) 745 (F)
		7.7 (22Q) 10.7 (22F)	11.9 (Q) 11.8 (F)		700				760 (Q) 785 (F)
100	100	6.7 (14Q) 8.1 (14F)	9.2	962 - 970	570	900	635	280	700
		6.9 (14Q) 8.3 (14F)	9.3		615			295	710
		6.7 (20Q) 9.7 (20F)	10.7	960 - 975	595		720	280	700
		6.9 (20Q) 9.8 (20F)	10.9		645			295	710
100	100	8.6	8.6	962 - 970	560	905	852	245	710
		8.8	8.8	962 - 970	610				
105	115	8.2 / 8.6 (14Q) 10.3 / 10.7 (14F)	12.1 / 12.5	1010 - 1030	575	885 / 935	725 (Q) 725 (F)	245 (Q) 285 (F)	640/660 (Q) 665/715 (F)
		8.5 / 8.9 (14Q) 10.6 / 11.0(14F)	12.4 / 12.8		625				690/710 (Q) 665/715 (F)
		8.8 / 9.2 (14Q) 10.9/ 11.3 (14F)	12.7 / 13.1		675				740/760 (Q) 720/760 (F)
	115 / 125	8.4 / 8.6 (20Q) 11.5 / 11.7 (20F)	13.3 / 13.5		595		725 (Q) 780 (F)	245 (Q) 310 (F)	640/660 (Q) 665/715 (F)
		8.7 / 8.9 (20Q) 11.8 / 12.0 (20F)	13.6 / 13.8		645				690/710 (Q) 665/715 (F)
		9.0 / 9.2 (20Q) 12.1 / 12.3 (20F)	13.9 / 14.1		695				740/760 (Q) 720/760 (F)
		8.6 (24Q) 12.2 (24F)	14.0		935	725 (Q) 825 (F)	660 (Q) 715 (F)		
		8.9 (24Q) 12.5 (24F)	14.3				710 (Q) 715 (F)		
		9.2 (24Q) 12.8 (24F)	14.6				695		760 (Q) 760 (F)
105	115	8.0 / 8.4 (14Q) 10.1 / 10.5 (14F)	11.8 / 12.2	1010 - 1030	575	885 / 935	725 (Q) 725 (F)	245 (Q) 285 (F)	640/660 (Q) 665/715 (F)
		8.3 / 8.7 (14Q) 10.4 / 10.8 (14F)	12.1 / 12.5		625				690/710 (Q) 665/715 (F)
		8.6 / 9.0 (14Q) 10.7/ 11.1 (14F)	12.4 / 12.8		675				740/760 (Q) 720/760 (F)
	115 / 125	8.2 / 8.4 (20Q) 11.3 / 11.5 (20F)	13.0 / 13.2		595		725 (Q) 780 (F)	245 (Q) 310 (F)	640/660 (Q) 665/715 (F)
		8.5 / 8.7 (20Q) 11.6 / 11.8 (20F)	13.3 / 13.5		645				690/710 (Q) 665/715 (F)
		8.8 / 9.0 (20Q) 11.9 / 12.1 (20F)	13.6 / 13.8		695				740/760 (Q) 720/760 (F)
		8.4 (24Q) 12.0 (24F)	13.7		935	725 (Q) 825 (F)	660 (Q) 715 (F)		
		8.7 (24Q) 12.3 (24F)	14.0				710 (Q) 715 (F)		
		9.0 (24Q) 12.6 (24F)	14.3				695		760 (Q) 760 (F)
97	115	9.1 (14Q) 11.3 (14F)	12.9	1010	575	905	735	310	705
		9.3 (14Q) 11.5 (14F)	13.1		625				
		9.5 (14Q) 11.7 (14F)	13.3		675				
		9.1 (24Q) 13.0 (24F)	14.6	1000	580	910	805	320	710
		9.4 (24Q) 13.2 (24F)	14.8		630				730
		9.6 (24Q) 13.4 (24F)	15.0		680				750
95		10.2(14Q) 12.7 (14F)	14.5	1090	630	975	845	350	755
		10.3 (14Q) 13.0 (14F)	14.8		680				780
		9.7 (24Q) 13.5 (24F)	15.3		655		875		775
		9.9 (24Q) 13.7 (24F)	15.5		705				800
95	115	8.1 (14Q) 10.3 (14F)	11.9	1010	575	905	735	310	705
		8.3 (14Q) 10.5 (14F)	12.1		625				
		8.5 (14Q) 10.7 (14F)	12.3		675				
		8.1 (24Q) 12.0 (24F)	13.6	1000	580	910	805	320	710
		8.4 (24Q) 12.2 (24F)	13.8		630				730
		8.6 (24Q) 12.4 (24F)	14.0		680				750

Series	Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)
S-Ergo	S-Ergo 115	440	470	7 / 14	515	425	230	360 - 435	14	8
		500		7 / 24						9
		440								
		500								
	S-Ergo 105	400	430	7 / 14	500	430	230	360 - 435	13	7
		450		7 / 24						
		400								
		450								
	S-Ergo 305	407	432	24"	495-515	400 (15.7 in)	220-270 (8.7 in-10.6 in)	415 - 490	-	-
		458								
Classic	KM-2500	420	400	6 / 14	450	420	225	310 - 410	11	9
		470								
	KM-2500L	420	400	6 / 20	450	420	225	310 - 410	11	9
		470								
	KM-1500	410	415	8 / 24	485	380	225	295 - 415	7	7
		460								
Flexx	Flexx	350	410 - 460	Caster: 7/ 8 Rear wheel: 14 / 24	390 - 530	430	200 - 280	330 - 430 Rear wheel 24 330 - 430 Rear wheel 14	7 - 27	5
		375								
		400								
		425								
		450								
		475								
		500								
	Flexx S (Self-propel)	350	415 - 465	Caster: 7/ 8 Rear wheel: 24	390 - 530	430	225 - 245	315 - 415	7 - 27	5
		375								
		400								
		425								
		450								
		475								
		500								
	Flexx S (Transport)	350	415 - 465	Caster: 7/ 8 Rear wheel: 14	390 - 530	430	225 - 245	315 - 415	7 - 27	5
		375								
		400								
		425								
		450								
		475								
		500								
	Flexx Tall	395 / 420 / 445 / 470	510 / 535	7/8 x 24	460 - 530	430 / 510	205 - 285	330 - 430	11	5 - 8
		395 / 420 / 445 / 470		7/8 x 14						5 - 7
	Flexx HD	500	390 - 490	7/8 x 24	460 - 530	430 / 510	215 - 295	370 - 475	12	6
		550								
		600								
		500		7/8 x 14			210 - 290	360 - 470	13	7
		550								
		600								

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)
95	115	9.3 (14Q) 11.7 (14F)	13.5	1090	615	975	820	330	740
		9.5 (14Q) 11.9 (14F)	13.7		665				780
		9.1 (24Q) 13.0 (24F)	14.8		640		855	315	760
		9.3 (24Q) 13.2 (24F)	15.0		690				800
98	115	11.5	12.1	1010	585	900	910	325	735
		11.7	12.3		635				
		9.4 (24Q) 13.2 (24F)	13.8	1005	590	905	980	335	740
		9.6 (24Q) 13.4 (24F)	14.0		640				
-	115	-	14.44	1080	595	970	728	320	800
			14.65		645		810	330	
101	100	10.7	11.2	885 - 915	550	855	740	250	680
		10.8	11.3		600				
101	100	11.7	12.2	900 - 930	575	855	790	285	705
		11.8	12.3		625				
108	100	13.8	14.1	1010 - 1060	720	865	915	350	770
		13.9	14.2		770				
99	130	9.8 - 10.9 Rear wheel 24 12.5 - 13.4 Rear wheel 14F 9.8 - 10.7 Rear wheel 14Q	15.9 - 16.7 Rear wheel 24 14.9 - 16.0 Rear wheel 14	1030 - 1105 Rear wheel 24 1005 - 1075 Rear wheel 14	560 / 525	915 - 940 Rear wheel 24 910 - 940 Rear wheel 14	815 - 865 Rear wheel 24 815 - 1040 Rear wheel 14	290	610 - 740
					585 / 550			290	
					610 / 575			305	
					635 / 600			305	
					660 / 625			320	
					685 / 650			320	
					710 / 675			335	
98	130	10.0 ~ 10.7	15.1 - 15.9	1025 - 1105	555	900 - 925	815 - 865	290	610 - 730
					580			305	
					605				
					630			335	
					655				
					680				
					705				
98	130	12.7 - 13.4 (14F) 10.0 - 10.7 (14Q)	14.1 - 14.9	1025 - 1105	520	900 - 925	815 - 1040	290	610 - 730
					545			305	
					570				
					595			335	
					620				
					645				
					670				
98	130	10.5 - 11.1	17.5 - 18.1	1130 - 1165	605 / 630 / 655 / 680	1000 - 1040	955 - 1115	340	700 - 950
			16.7 - 17.3		645 / 670 / 695 / 720	975 - 1000	855 - 1115		
98	170	11.3 - 12.1	18.1 - 19.5	1100 - 1220	720	930 - 960	1050 - 1130	320 - 330	740 - 835
					770				
					820				
		14.1 - 14.9	17.4 - 18.8	1080 - 1200	680	920 - 950			
					730				
					780				

Series	Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)
Flexx	Flexx Adapt	255 - 355	230 - 380	6 / 16	480 - 530	350 - 500	10 - 230	160 - 360	0 - 40	0 - 45
		255 - 355		8 / 16						
	Flexx Junior	255 - 355	285 - 365	5/6 x 14	410	380	155 - 235	150 - 250	11 - 31	5
		255 - 355		5/6 x 20	400		180 - 235		13 - 33	
		255 - 355		5/6 x 22	415		220 - 235		9 - 29	2
Positioning	SME	400	400	6 x 24	490	415	160 - 240	350 - 410	10	1
		440	440							
		480	480							
	VIP 2	400	415	6/7 x 14	470 / 490	510 / 560	220 - 280	320 - 420	13 - 43	1 - 36
		400	465							
		450	415							
		450	465							
	VIP 2	400	415	6 / 7 x 20	470 / 490	510 / 560	220 - 280	320 - 420	13 - 43	1 - 36
		400	465							
		450	415							
		450	465							
	VIP 515	405	415	7 x 14	490	460	240 - 320	405 - 480	13	1 - 37
		455			485					
		405		7 x 20	490					
		455			485					
	MVP 502	420	420	7 x 14	520	490	260 - 335	410 - 485	13 - 64	7
		420								
		450								
		450								
		475	470		525					
		500								
		420	420	7 x 22	520					
		420								
		450								
		450								
		475	470							
		500								
	Ergo Recline	415	420	8 x 14	500	495	220 - 300	365 - 440	8 - 80	9
		465								
		415		8 x 20						
		465								
	KM-5000	425	420	8 x 14	500	450	235	450 - 550	7 - 76	7
		460	470				230		375 - 425	
		425	420							
		460	470							

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)
100	75	18.6	29.8	1070 - 1385	570	1060 - 1190	770	570	530
94	75	9.7 - 9.8	14.3 - 14.4	820 - 835	450 - 550	800 - 850	810 - 870	325	725
		9.4 - 9.5	14.9 - 15.0	830 - 900	485 - 590	780 - 830		355	
		9.4 - 9.5	15.2 - 15.3	840 - 890	485 - 590	825 - 875		355	
92	120	18.8	27.4	1005 - 1060	590	950	1015	535	600
		21.0	29.8		630			580	
		23.2	32.0		670			625	
100	115	17.8	24.4	1070	630	1380	775	355	680
		18.0	24.7	1120			825		
		18.1	24.9	1070	680		775		
		18.3	25.2	1120			825		
100	115	15.4(20Q) 19.0(20F)	25.4	1070	660	1380	840	380	715
		15.6(20Q) 19.2(20F)	25.6	1120			890		
		15.6(20Q) 19.2(20F)	25.6	1070	710		840		
		15.8(20Q) 19.4(20F)	25.8	1120			890		
114	115	15.1	18.9	1070	635	955 / 1325 (with headrest)	725	350	760
		15.3	19.3		680		730		
		16.2	20.0		660	955 / 1330 (with headrest)	790	375	
		16.4	20.3		710				
100 - 160	115	14.5	20.3	1190	550	840 / 1350 (with headrest)	785	305	770
		14.6	20.4		580				
		14.7	20.5		610				
		14.8	20.8		635				
		15.0	20.9	1225	635		825		
		15.4	21.0		700				
		16.8	21.9	1190	580	840 / 1350 (with headrest)	870	360	770
		16.8	21.9		610				
		16.9	22.0		640				
		17.0	22.3		665				
		17.2	22.5	1225	665		905		
		17.6	22.9		730				
115 - 192	115	13.3	21.1	1190 - 1815	590	700 ~ 1300	790	340	940
		13.4	21.2		640				
		14.2	21.8		640		850	400	
		14.3	21.9		690				
114 - 180	115	13.2	18.7	1215	660	1295	745	330	935
13.5		19.0	1265	795					
110 - 117		15.7	21.3	1240			910	400	930
		16.0	21.6	1290			960		

Series	Product Name	Nominal seat width (mm/in)	Nominal seat depth (mm/in)	Front caster / Rear wheel size (in)	Seat surface height at front edge (mm)	Backrest height (mm)	Armrest to seat distance (mm)	Footrest to seat distance (mm)	Backrest angle (X°)	Seat plane angle (X°)								
Active & Specialty	KM-BT10	560/610/660/710/760	460/510/560/610	7 / 8 x 24Q	390 - 500	410 - 510	250	320 - 370	(-5) - 15	-								
	KM-AT60	380 - 455	355 - 505	5 x 24Q	460 - 580	310 - 480	-	-	-	0/3/5								
	Ergo Live	340 / 14	390 / 15	5 x 24	465	360 - 455	225 / 250	375 - 435	7 / 12 / 17	5								
		340 / 14	415 / 16															
		390 / 16	415 / 16															
		390 / 16	465 / 18															
		440 / 18	415 / 16															
440 / 18	465 / 18																	
SOMA	SOMA 215	425	405	6 x 14	430	400	220	345 - 400	11	8								
		425		6 x 20						6								
		460	430	7 x 14	480	430		430 - 474	10	5								
		460		7 x 20														
		460		7 x 22Q														
	SOMA 100	400	440	7 x 22	480	390	220	395 - 445	11	5								
		450																
	SOMA 105	400	440	7 x 22	480	390	220	395 - 445	11	5								
		450																
	SOMA 103	420	430	7 x 16	500	400	235	395 - 420	13	5								
				7 x 22														
	AGILE		16	16	7 x 14	480	430 / 510	220	380 - 480	14	8							
			18															
			16	18								7 x F22	480	430 / 510	220	380 - 480	14	8
			18															
			20															
			16	16	8 x 14	505	430 / 510	220	380 - 480	14	8							
			18															
			16	18								8 x Q24	505	430 / 510	220	380 - 480	14	8
			18															
			20															
			16	16	8 x Q24	505	430 / 510	220	380 - 480	14	8							
			18															
			16	18								8 x Q24	505	430 / 510	220	380 - 480	14	8
			18															
20																		

Leg to seat angle (X°)	Maximum user weight (kg)	Mass of the heaviest part (kg)	Total mass (kg)	Overall length with legrest (mm)	Overall width (mm)	Overall height (mm)	Folded length (mm)	Folded width (mm)	Folded height (mm)
-	250	22.9	31	Min.1180	Min.780	Min.800	Min.900	Min.440	Min.800
				Max.1230	Max.980	Max.1010	Max.950	Max.440	Max.1010
-	100	Min. 8	Min. 10.5	880 - 1030	620 - 700	770 - 1060	-	-	-
97	125	9.0	13.1	895 - 910	580	810 - 885	855	290	550
		9.1	13.2	920 - 935			870		
		9.1	13.3	920 - 935	630		870	310	585
		9.2	13.5	970 - 985			900		
		9.2	13.5	920 - 935	680		870	330	620
		9.3	13.7	970 - 985			900		
95	100	10.2	10.2	940	605	855	630	300	690
		11.2	11.2	950			715	310	
110		10.5	10.9	1010	640	940	680	255	760
		11.7	12	1020			775	280	
		8.5	12.4	1045			660	785	
96	100	11.1	11.3	975	610	875	955	285	875
		11.4	11.6		665				
96	100	11.2	11.4	975	610	875	955	285	760
		11.5	11.7		665				
100	100	10.4	10.4	945	585	875	945	280	865
		11.6	11.6	980	615	865	960	295	
116	125	11.1 - 11.6 11.4 - 11.9 11.6 - 12.1	12.3 - 12.8 12.6 - 13.1 12.8 - 13.3	1025 - 1035	600	915 / 995	710 / 780	280	735 / 915
					650				
				1085 - 1095	600		760 / 830		
					650				
					700				
116	125	11.7 - 12.2 12.0 - 12.5 12.2 - 12.7	12.9 - 13.4 13.2 - 13.7 13.4 - 13.9	1055 - 1075	595	920 / 1000	710 / 770	265	765 / 920
					645				
				1105 - 1125	595		760 / 820		
					645				
					695				
116	125	11.0 - 11.5 11.3 - 11.8 11.5 - 12.0	12.2 - 12.7 12.5 - 13.0 12.7 - 13.2	1025 - 1035	600	935 / 1015	750 / 790	280	760 / 930
					650				
				1085 - 1095	600		800 / 840		
					650				
					700				
116	125	8.6 - 9.1 8.9 - 9.4 9.1 - 9.6	13.9 - 14.4 14.2 - 14.7 14.4 - 14.9	1065 - 1070	605	940 / 1020	850 / 1020	280	750
					655				
				1115 - 1120	605		900 / 1070		
					655				
					705				