

Product Name	Nominal seat width (mm)	Nominal seat depth (mm)	Effective seat width (mm)	Effective seat depth (mm)	Backrest height w/cushion (mm)	Seat height w/o cushion (mm)	Seat height w/cushion (mm)	Power seat lift (mm)	Power tilt, seat plane (X°)	Backrest recline, seat to back angle (X°)	Legrest angle, leg to seat (X°)
EVO Altus	420/470/520	400 - 600	420 - 520	400 - 600	510 - 670	415 - 865	490 - 940	400	0 - 50	85 - 180	85 - 180
EVO Lectus LR	420/470/520	400 - 600	420 - 520	400 - 600	510 - 670	415 - 865	480 - 930	400	3 - 48	85 - 180	85 - 180
Morgan M Lectus	420/470/520	425 - 600	420 - 520	425 - 600	570 - 600	485 - 785	590 - 890	300	3 - 48	80 - 177	98 - 173
Morgan M Versus w/ Lectus cushion	420/470/520	450/525	340 - 665	150 - 550	525 - 765	440 - 740	540 - 840	300	3 - 48	90 - 135	105 - 185
Morgan M Sling	430/480	400 - 475 475 - 500	430 - 520 450 - 520	400 - 475 475 - 500	465/545	420	500	-	2 - 20	90 - 122	105
Leon F Versus w/ Lectus cushion	420/470/520	450/525	340 - 665	150 - 550	525 - 765	440 - 740	540 - 840	300	2 - 47	90 - 135	105 - 185
Leon R Versus w/ Lectus cushion	420/470/520	450/525	340 - 665	150 - 550	525 - 765	460 - 760	560 - 860	300	4 - 49	90 - 135	105 - 185
Leon R Sling	430/480	400 - 475 475 - 500	430 - 520 450 - 520	400 - 475 475 - 500	465/545	460	540	-	3 - 21	90 - 122	109
Leon F Kameleon Size 1	250 - 300	255 - 395	250 - 300	270 - 410	390 - 450	463 - 763	503 - 803	300	5 - 45	85 - 130	90 - 140
Leon F Kameleon Size 2	350 - 400	310 - 450	350 - 400	330 - 470	455 - 550	463 - 763	503 - 803	300	5 - 45	85 - 130	90 - 140
Blazer Sling	430/480	375 - 475 425 - 475	420/470	375 - 475 425 - 475	460/540	445	540	-	8 - 28	92/97/102	105
Blazer Captain	450/515	460/480	470/500	420 - 500 490 - 570	535	-	565/553	-	-	0 - 134 0 - 131	100/105
Saber Sling	430	400 - 475	450 - 535	400 - 475	455	525	515	-	3 - 36	96 - 133	117
	480	475 - 500	485 - 535	475 - 500	520	600	515	-	3 - 36	96 - 133	117
Saber Captain	450	460	470 - 525	460	520	-	520 550 (T)	-	4 - 38	0 - 135 (Manual) 91 - 121 (power)	104
	515	480	510 - 525	480	540	-	520 550 (T)	-	8 - 42	0 - 135 (Manual) 91 - 121 (power)	104
Ergo Stand	400/440/480	460/490/520	420/455/490	500/530/560	455	-	495	-	-	75 - 115	92/94/96
eFlexx-14	360/410/460	415/440/465	405/455/500	560/585/610	410	-	485	-	-	90 - 110	99
eFlexx-20	360/410/460	415/440/465	390/455/510	545/585/610	410	-	485	-	-	90 - 110	99
Falcon	410/460/515	430/480/480	430/480/530	470/515/515	450-525	-	525	-	-	95	110
Ergo Traveler	445	460	450	445	500	-	540	-	-	88 - 102	90
KS-343.2S	425 (17")	420 (17")	425	420	490	-	592 - 642 (10"wheel) 600 - 650 (11"wheel)	-	-	101.5	90
KS-747.2	430 (18")	435 (18")	430	435	510	-	610 - 635	-	-	90 - 130	90
KS-737.2	430 (18")	435 (18")	430	435	510	-	600 - 625	-	-	90 - 130	90

* Drive range is estimated based on the following conditions: 20°C - 35°C, 100 kg occupant with VR2 system / 136 kg occupant with LiNX system, brand new fully charged batteries and driving on the flat road.

F/M/R Wheel (in)	Overall length with legrest (mm)	Overall width (mm)	Max. User weight (kg)	Drive range (km)	Max. speed (km)	Obstacle climbing (mm)	Min. turning radius w/o user (mm)	Ability to climb rated slope (X°)	Battery (Ah)	Electrical system
14/-/9	1200	610 - 660	136	54	12	70	750	12	80*2	R-net
14/-/9	1200	610 - 660	136	54	12	70	750	12	80*2	R-net
8/14/6	1240 - 2120	620 - 700	136	33	12	75	835	12	80*2	R-net
8/14/6	1195 - 1660	600 - 825	136	39.7	12	75	805	12	80*2	VR2/R-net
8/14/6	1090 - 1120	630 - 690	136	27 (50Ah) 49 (80Ah)	12	75	700	12	50*2 / 80*2	VR2/R-net
14/-/9	1130 - 1400	620 - 825	136	39	12	65	660	12	80*2	VR2/R-net
9/-/14	1200 - 1550	620 - 825	136	35.6	12	50	950	12	80*2	VR2/R-net
9/-/14	1140 - 1180	630 - 690	136	44	12	50	895	12	80*2	VR2/R-net
14/-/9	1018	625	57	25.6 (50 Ah) 46.2 (80 Ah)	12	50	700	12	50*2 / 80*2	VR2/R-net
	1120									
9/-/14	1090 - 1160	580/600	136	27 (VR2)* 30 (Linx)*	10	50	825	12	50*2	VR2/Linx
9/-/14	1060 - 1130	600/625	136	27 (VR2)* 30 (Linx)*	10	50	825	12	50*2	VR2/Linx
9/-/16	1130 - 1185	626 - 680	136	26.5 (50 Ah) 42.4 (80 Ah)	10 / 13	65	905	12	50*2 / 80*2	VR2
9/-/16	1130 - 1185	626 - 680	136		10 / 13	65	905	12		
9/-/16	1135 - 1545	626 - 680	136	26.5 (50Ah) 42.4 (80Ah)	10 / 13	65	905	12	50*2 / 80*2	VR2
9/-/16	1135 - 1650	665 - 680	136		10 / 13	65	905	12	50*2	VR2
12.5 /-/6	1000 - 1120/ 1005 - 1165/ 1005 - 1205	640	135	25.9	9.3	60	745	12	50*2	VR2/Linx
8/-/14	1015 - 1065	560/615/ 665	100	10 (11.5 AH) 20 (26.1 AH)	8	35	795/825/ 905	10	11.5*1 26.1*1	Linx
8/-/20	1125 - 1175	590/645/ 690	100	10 (11.5 AH) 20 (26.1 AH)	8	35	795/840/ 865	10	11.5*1 26.1*1	Linx
8/-/12	985 - 1115	590/640/ 690	115/135 (pneumatic tires)	32	6	50	650/825/ 1000	12	50*2	VR2
8/-/12	1030 - 1075	620	120	14.9 (22 Ah) 18.7 (36 Ah)	6	50	755	6	22*2 / 36*2	VR2
10/-/10 11/-/11	1180 - 1280	600-680	136	25 - 28 (36 Ah) 32 - 37 (36 Ah)	9 (10" wheel) 10 (11" wheel)	75	1600	10	36*2 / 50*2	S-Drive
13/-/13	1330 - 1800	670	135	42	15	100	1575	10	50*2	S-Drive
10/-/13	1230 - 1650	615	135	42	15	75	1090	10	50*2	S-Drive

- The actual product specifications may differ from the data in the table by ± 10 mm ± 0.5 kg.
- The actual product specifications may vary according to different configurations.
- Karma reserves the right to modify information herein without further notice.