



A dynamic partner that goes everywhere with you

A sleek design that breaks the stereotypical look of traditional mobility aids. iBuddy is a great companion for the outdoors, designed to fit elegantly into the home environment, with unobstructed access indoors and outdoors.

Excellent Outdoor Performance


iBuddy comes equipped with four-wheel independent suspension for a smoother ride as well as powerful motors for superior handling over rough terrain. Up to 75 mm of obstacle climbing to easily overcome thresholds and bumpy roads. In addition, LED lighting provides better visibility while you travel.

Surroundings Made Accessible

iBuddy has a power sliding forward/backward seat function, which help user to sit down with ease and access the desk or dinner table, making socialization more delightful.

Extend Your Daily Routine

With a range of up to 50 km, you can shop for groceries, go to school, and make friends without worrying about your wheelchair.



Equip your iBuddy with an iSeat smart cushion, and maintain good sitting habits starting today!



karma
A Better Wheelchair, A Better Fit



KARMA MEDICAL PRODUCTS CO., LTD.
Email: globalsales@karma.com.tw
www.KarmaMedical.com